



**Blue plaster material**

Can be made from washproof PE, washproof and elastic PU, 2nd skin film or elastic textile material. Blue color for visual detection

**Strong hypo-allergenic adhesive**

Strong adhesive ensures to stay in place in difficult work conditions, even in contact with water

**Aluminium strip**

For detection by metal detectors in the food industry

**X Ray thread**

Visible for X ray detectors

**Absorbent wound pad**

**Non adherent wound contact layer**

**Protective foil**

To protect against dirt and dust

**Wrapper**

