

SPECIAL BLUE COLOUR PRODUCTS FOR HOTEL, CATERING AND INDUSTRIAL KITCHENS



Non Woven Cleaning and Sanitizing Towel

Cleaning & Sanitizing ready to use towels are impregnated with disinfectant for efficient & deep cleaning of hands & small surfaces. They are alcohol free hence no degradation of material. Acc to HACCP, BRC & IFS.

Content : 100 Towels per Bucket
Packaging : 12 buckets/carton



Hydrogel Burn Plasters for Minor Burns & Abrasions

The sterile Hydrogel contains 60% water which has a slight cooling effect and reduces redness and immediate pain. The PU Film allows wound to breathe & obstruct bacteria. Aluminum Strip underneath wound pad for detection by metal detectors in food industry.

Packaging : 10 Pcs/box.



Self Adhesive Plaster

Plaster for small injuries & distortion. Made from non woven polyester material which is breathable & does not stick to wound. Blue in color for visual detection. Acc to HAACP & BRC.

CONMED DEVICES PVT. LTD.

Corporate Office

42, Todarmal Road, Bengali Market, Near Mandi House, New Delhi 110 001, India

P: +91 11 43512874 | M: +91 98102 44134, 84477 51030 | E: ab@conmed.in, db@conmed.in | W: www.conmed.in